



FKTN Antibody, Biotin conjugated

Product Code	CSB-PA008709LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75072
Immunogen	Recombinant Human Fukutin protein (87-276AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Fukutin (EC 2.-.-.-) (Fukuyama-type congenital muscular dystrophy protein), FKTN, FCMD
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	FKTN
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.