



RACK1 Antibody, FITC conjugated

Product Code	CSB-PA009605LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P63244
Immunogen	Recombinant Human Receptor of activated protein C kinase 1 protein (93-317AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Receptor of activated protein C kinase 1 (Cell proliferation-inducing gene 21 protein) (Guanine nucleotide-binding protein subunit beta-2-like 1) (Guanine nucleotide-binding protein subunit beta-like protein 12.3) (Human lung cancer oncogene 7 protein) (HLC-7) (Receptor for activated C kinase) (Small ribosomal subunit protein RACK1) [Cleaved into: Receptor of activated protein C kinase 1, N-terminally processed (Guanine nucleotide-binding protein subunit beta-2-like 1, N-terminally processed)], RACK1, GNB2L1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	RACK1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.