



HCRTR1 Antibody, HRP conjugated

Product Code	CSB-PA010231HB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43613
Immunogen	Recombinant Human Orexin receptor type 1 protein (1-46AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Orexin receptor type 1 (Ox-1-R) (Ox1-R) (Ox1R) (Hypocretin receptor type 1), HCRTR1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	HCRTR1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.