



HHEX Antibody, Biotin conjugated

Product Code	CSB-PA010335OD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q03014
Immunogen	Peptide sequence from Human Hematopoietically-expressed homeobox protein HHEX protein (187-204AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Hematopoietically-expressed homeobox protein HHEX (Homeobox protein HEX) (Homeobox protein PRH), HHEX, HEX PRH PRHX
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	HHEX
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.