



IFNA5 Antibody, Biotin conjugated

Product Code	CSB-PA011041ND01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01569
Immunogen	Peptide sequence from Human Interferon alpha-5 protein (52-73AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Interferon alpha-5 (IFN-alpha-5) (Interferon alpha-61) (Interferon alpha-G) (LeIF G), IFNA5
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	IFNA5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.