



IFNA5 Antibody, Biotin conjugated

| | |
|---------------------|---|
| Product Code | CSB-PA011041ND01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P01569 |
| Immunogen | Peptide sequence from Human Interferon alpha-5 protein (52-73AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Biotin |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | Antigen Affinity Purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Interferon alpha-5 (IFN-alpha-5) (Interferon alpha-61) (Interferon alpha-G) (LeIF G), IFNA5 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Immunology |
| Target Names | IFNA5 |