



ITGB1BP1 Antibody, FITC conjugated

Product Code	CSB-PA011881LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O14713
Immunogen	Recombinant Human Integrin beta-1-binding protein 1 protein (1-200AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Integrin beta-1-binding protein 1 (Integrin cytoplasmic domain-associated protein 1) (ICAP-1), ITGB1BP1, ICAP1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ITGB1BP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.