



MAP3K2 Antibody, Biotin conjugated

Product Code	CSB-PA013424LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y2U5
Immunogen	Recombinant Human Mitogen-activated protein kinase kinase kinase 2 protein (357-619AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mitogen-activated protein kinase kinase kinase 2 (EC 2.7.11.25) (MAPK/ERK kinase kinase 2) (MEK kinase 2) (MEKK 2), MAP3K2, MAPKKK2 MEKK2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MAP3K2