



MAP3K2 Antibody, Biotin conjugated

Product Code	CSB-PA013424LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y2U5
Immunogen	Recombinant Human Mitogen-activated protein kinase kinase kinase 2 protein (357-619AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mitogen-activated protein kinase kinase kinase 2 (EC 2.7.11.25) (MAPK/ERK kinase kinase 2) (MEK kinase 2) (MEKK 2), MAP3K2, MAPKKK2 MEKK2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MAP3K2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.