



MEF2A Antibody, FITC conjugated

Product Code	CSB-PA013670LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q02078
Immunogen	Recombinant Human Myocyte-specific enhancer factor 2A protein (72-226AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Myocyte-specific enhancer factor 2A (Serum response factor-like protein 1), MEF2A, MEF2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MEF2A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.