



PRAMEF18 Antibody, HRP conjugated

Product Code	CSB-PA018613LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5VWM3
Immunogen	Recombinant Human PRAME family member 18 protein (1-257AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	PRAME family member 18, PRAMEF18
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	PRAMEF18
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.