



PSMB8 Antibody, FITC conjugated

Product Code	CSB-PA018886LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P28062
Immunogen	Recombinant Human Proteasome subunit beta type-8 protein (176-276AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Proteasome subunit beta type-8 (EC 3.4.25.1) (Low molecular mass protein 7) (Macropain subunit C13) (Multicatalytic endopeptidase complex subunit C13) (Proteasome component C13) (Proteasome subunit beta-5i) (Really interesting new gene 10 protein), PSMB8, LMP7 PSMB5i RING10 Y2
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	PSMB8
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.