



# RHOB Antibody, Biotin conjugated

Product Code	CSB-PA019683LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P62745
Immunogen	Recombinant Human Rho-related GTP-binding protein RhoB protein (1-196AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Rho-related GTP-binding protein RhoB (Rho cDNA clone 6) (h6), RHOB, ARH6, ARHB
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	RHOB