



ROS1 Antibody, Biotin conjugated

Product Code	CSB-PA020072LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P08922
Immunogen	Recombinant Human Proto-oncogene tyrosine-protein kinase ROS protein (79-248AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Proto-oncogene tyrosine-protein kinase ROS (EC 2.7.10.1) (Proto-oncogene c-Ros) (Proto-oncogene c-Ros-1) (Receptor tyrosine kinase c-ros oncogene 1) (c-Ros receptor tyrosine kinase), ROS1, MCF3 ROS
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ROS1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.