



RTTN Antibody, FITC conjugated

Product Code	CSB-PA020581LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86VV8
Immunogen	Recombinant Human Rotatin protein (1509-1603AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Rotatin, RTTN
Immunogen Species	Homo sapiens (Human)
Research Area	Developmental Biology
Target Names	RTTN
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.