



SBK2 Antibody, Biotin conjugated

Product Code	CSB-PA020731LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0C263
Immunogen	Recombinant Human Serine/threonine-protein kinase SBK2 protein (62-330AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serine/threonine-protein kinase SBK2 (EC 2.7.11.1) (SH3-binding domain kinase family member 2) (Sugen kinase 69) (SgK069), SBK2, SGK069
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SBK2