



SBK2 Antibody, Biotin conjugated

Product Code	CSB-PA020731LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0C263
Immunogen	Recombinant Human Serine/threonine-protein kinase SBK2 protein (62-330AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serine/threonine-protein kinase SBK2 (EC 2.7.11.1) (SH3-binding domain kinase family member 2) (Sugen kinase 69) (SgK069), SBK2, SGK069
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SBK2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.