



SOX10 Antibody, Biotin conjugated

Product Code	CSB-PA022418OD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P56693
Immunogen	Peptide sequence from Human Transcription factor SOX-10 protein (49-68AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transcription factor SOX-10, SOX10
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	SOX10