



TBKBP1 Antibody, Biotin conjugated

Product Code	CSB-PA023233LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A7MCY6
Immunogen	Recombinant Human TANK-binding kinase 1-binding protein 1 protein (327-427AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	TANK-binding kinase 1-binding protein 1 (TBK1-binding protein 1), TBKBP1, KIAA0775 SINTBAD
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	TBKBP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.