



TRO Antibody, Biotin conjugated

Product Code	CSB-PA025069LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q12816
Immunogen	Recombinant Human Trophinin protein (1-164AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Trophinin (MAGE-D3 antigen), TRO, KIAA1114 MAGED3
Immunogen Species	Homo sapiens (Human)
Research Area	Cancer
Target Names	TRO
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.