



TYRP1 Antibody, HRP conjugated

Product Code	CSB-PA025399LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P17643
Immunogen	Recombinant Human 5,6-dihydroxyindole-2-carboxylic acid oxidase protein (91-186AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	5,6-dihydroxyindole-2-carboxylic acid oxidase (DHICA oxidase) (EC 1.14.18.-) (Catalase B) (Glycoprotein 75) (Melanoma antigen gp75) (Tyrosinase-related protein 1) (TRP) (TRP-1) (TRP1), TYRP1, CAS2 TYRP TYRRP
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	TYRP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.