




USP46 Antibody

Product Code	CSB-PA025738LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P62068
Immunogen	Recombinant Human Ubiquitin carboxyl-terminal hydrolase 46 protein (66-166AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IF; Recommended dilution: IF:1:200-1:500
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ubiquitin carboxyl-terminal hydrolase 46 (EC 3.4.19.12) (Deubiquitinating enzyme 46) (Ubiquitin thioesterase 46) (Ubiquitin-specific-processing protease 46), USP46
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	USP46

Image



Immunofluorescence staining of A549 cells with CSB-PA025738LA01HU at 1:267, counter-stained with DAPI. The cells were fixed in 4% formaldehyde, permeabilized using 0.2% Triton X-100 and blocked in 10% normal Goat Serum. The cells were then incubated with the antibody overnight at 4°C. The secondary antibody was Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).