



AOX1 Antibody, HRP conjugated

Product Code	CSB-PA24692B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O65000
Immunogen	Peptide sequence from Chlamydomonas reinhardtii Ubiquinol oxidase protein (202-217AA)
Raised In	Rabbit
Species Reactivity	Chlamydomonas reinhardtii
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ubiquinol oxidase, AOX1
Immunogen Species	Chlamydomonas reinhardtii (Chlamydomonas smithii)
Research Area	Others
Target Names	AOX1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.