



SKOR1 Antibody, HRP conjugated

Product Code	CSB-PA308064HB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P84550
Immunogen	Recombinant Human SKI family transcriptional corepressor 1 protein (744-965AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	SKI family transcriptional corepressor 1 (Functional Smad-suppressing element on chromosome 15) (Fussel-15) (LBX1 corepressor 1) (Ladybird homeobox corepressor 1), SKOR1, CORL1 FUSSEL15 LBXCOR1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SKOR1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.