



# SKOR1 Antibody, Biotin conjugated

Product Code	CSB-PA308064HD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P84550
Immunogen	Recombinant Human SKI family transcriptional corepressor 1 protein (744-965AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	SKI family transcriptional corepressor 1 (Functional Smad-suppressing element on chromosome 15) (Fussel-15) (LBX1 corepressor 1) (Ladybird homeobox corepressor 1), SKOR1, CORL1 FUSSEL15 LBXCOR1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SKOR1