



# At1g32060 Antibody, FITC conjugated

<b>Product Code</b>	CSB-PA326237OC01DOA
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P25697
<b>Immunogen</b>	Peptide sequence from Arabidopsis thaliana Phosphoribulokinase, chloroplastic protein (223-235AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Arabidopsis thaliana
<b>Form</b>	Liquid
<b>Conjugate</b>	FITC
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Phosphoribulokinase, chloroplastic (PRK) (PRKase) (EC 2.7.1.19) (Phosphopentokinase)
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Research Area</b>	Others
<b>Target Names</b>	At1g32060
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.