



At1g32060 Antibody, Biotin conjugated

Product Code	CSB-PA326237OD01DOA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P25697
Immunogen	Peptide sequence from Arabidopsis thaliana Phosphoribulokinase, chloroplastic protein (223-235AA)
Raised In	Rabbit
Species Reactivity	Arabidopsis thaliana
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Phosphoribulokinase, chloroplastic (PRK) (PRKase) (EC 2.7.1.19) (Phosphopentokinase)
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Research Area	Others
Target Names	At1g32060
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.