



ppsA Antibody, HRP conjugated

Product Code	CSB-PA332921HB01ENV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P23538
Immunogen	Recombinant Escherichia coli Phosphoenolpyruvate synthase protein (2-792AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Phosphoenolpyruvate synthase (PEP synthase) (EC 2.7.9.2) (Pyruvate, water dikinase), ppsA, pps
Immunogen Species	Escherichia coli(strain K12)
Research Area	Others
Target Names	ppsA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.