



ppsA Antibody, Biotin conjugated

| | |
|---------------------|--|
| Product Code | CSB-PA332921HD01ENV |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P23538 |
| Immunogen | Recombinant Escherichia coli Phosphoenolpyruvate synthase protein (2-792AA) |
| Raised In | Rabbit |
| Species Reactivity | Escherichia coli |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Biotin |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Phosphoenolpyruvate synthase (PEP synthase) (EC 2.7.9.2) (Pyruvate, water dikinase), ppsA, pps |
| Immunogen Species | Escherichia coli(strain K12) |
| Research Area | Others |
| Target Names | ppsA |