



alaA Antibody, HRP conjugated

Product Code	CSB-PA359250LB01ENV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0A959
Immunogen	Recombinant Escherichia coli Glutamate-pyruvate aminotransferase AlaA protein (1-405AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glutamate-pyruvate aminotransferase AlaA (EC 2.6.1.2), alaA, yfbQ
Immunogen Species	Escherichia coli(strain K12)
Research Area	Others
Target Names	alaA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.