



JRK Antibody, Biotin conjugated

Product Code	CSB-PA525639ND01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75564
Immunogen	Peptide sequence from Human Jerky protein homolog protein (465-483AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Jerky protein homolog, JRK, JH8
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	JRK
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.