



RBL2 Antibody, Biotin conjugated

Product Code	CSB-PA600921LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q08999
Immunogen	Recombinant Human Retinoblastoma-like protein 2 protein (417-616AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Retinoblastoma-like protein 2 (130 kDa retinoblastoma-associated protein) (p130) (Retinoblastoma-related protein 2) (RBR-2) (pRb2), RBL2, RB2
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	RBL2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.