



ROCK1 Antibody, Biotin conjugated

Product Code	CSB-PA614415HD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13464
Immunogen	Recombinant Human Rho-associated protein kinase 1 protein (405-668AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Rho-associated protein kinase 1 (EC 2.7.11.1) (Renal carcinoma antigen NY-REN-35) (Rho-associated, coiled-coil-containing protein kinase 1) (Rho-associated, coiled-coil-containing protein kinase I) (ROCK-I) (p160 ROCK-1) (p160ROCK), ROCK1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ROCK1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.