



RABEP1 Antibody, HRP conjugated

Product Code	CSB-PA614886LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15276
Immunogen	Recombinant Human Rab GTPase-binding effector protein 1 protein (1-255AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Rab GTPase-binding effector protein 1 (Rabaptin-4) (Rabaptin-5) (Rabaptin-5alpha) (Renal carcinoma antigen NY-REN-17), RABEP1, RAB5EP RABPT5 RABPT5A
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	RABEP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.