



COL9A2 Antibody, FITC conjugated

Product Code	CSB-PA617900OC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14055
Immunogen	Peptide sequence from Human Collagen alpha-2(IX) chain protein (594-611AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Collagen alpha-2(IX) chain, COL9A2
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	COL9A2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.