



DUSP4 Antibody, FITC conjugated

Product Code	CSB-PA619752LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13115
Immunogen	Recombinant Human Dual specificity protein phosphatase 4 protein (109-198AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Dual specificity protein phosphatase 4 (EC 3.1.3.16) (EC 3.1.3.48) (Dual specificity protein phosphatase hVH2) (Mitogen-activated protein kinase phosphatase 2) (MAP kinase phosphatase 2) (MKP-2), DUSP4, MKP2 VH2
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	DUSP4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.