



LLGL1 Antibody, Biotin conjugated

Product Code	CSB-PA622993LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15334
Immunogen	Recombinant Human Lethal (2) giant larvae protein homolog 1 protein (953-1056AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Lethal(2) giant larvae protein homolog 1 (LLGL) (DLG4) (Hugl-1) (Human homolog to the D-lgl gene protein), LLGL1, DLG4 HUGL HUGL1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	LLGL1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.