



LLGL1 Antibody, Biotin conjugated

Product Code	CSB-PA622993LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15334
Immunogen	Recombinant Human Lethal (2) giant larvae protein homolog 1 protein (953-1056AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Lethal(2) giant larvae protein homolog 1 (LLGL) (DLG4) (Hugl-1) (Human homolog to the D-lgl gene protein), LLGL1, DLG4 HUGL HUGL1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	LLGL1