



NUGGC Antibody, Biotin conjugated

Product Code	CSB-PA718072LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q68CJ6
Immunogen	Recombinant Human Nuclear GTPase SLIP-GC protein (1-228AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nuclear GTPase SLIP-GC (EC 3.6.1.-) (Speckled-like pattern in the germinal center), NUGGC, C8orf80
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	NUGGC
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.