



NEK10 Antibody, HRP conjugated

Product Code	CSB-PA747816LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6ZWH5
Immunogen	Recombinant Human Serine/threonine-protein kinase Nek10 protein (1-104AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serine/threonine-protein kinase Nek10 (EC 2.7.11.1) (Never in mitosis A-related kinase 10) (NimA-related protein kinase 10), NEK10
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	NEK10