



TAOK1 Antibody, HRP conjugated

Product Code	CSB-PA752101LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7L7X3
Immunogen	Recombinant Human Serine/threonine-protein kinase TAO1 protein (400-659AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serine/threonine-protein kinase TAO1 (EC 2.7.11.1) (Kinase from chicken homolog B) (hKFC-B) (MARK Kinase) (MARKK) (Prostate-derived sterile 20-like kinase 2) (PSK-2) (PSK2) (Prostate-derived STE20-like kinase 2) (Thousand and one amino acid protein kinase 1) (TAOK1) (hTAOK1), TAOK1, KIAA1361 MAP3K16 MARKK
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TAOK1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.