



USP17L1 Antibody, HRP conjugated

Product Code	CSB-PA755915LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7RTZ2
Immunogen	Recombinant Human Ubiquitin carboxyl-terminal hydrolase 17-like protein 1 protein (19-235AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ubiquitin carboxyl-terminal hydrolase 17-like protein 1 (EC 3.4.19.12) (Deubiquitinating enzyme 17-like protein 1) (Ubiquitin thioesterase 17-like protein 1) (Ubiquitin-specific-processing protease 17-like protein 1), USP17L1, USP17L USP17L1P
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	USP17L1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.