



CD163 Antibody, HRP conjugated

Product Code	CSB-PA801238LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86VB7
Immunogen	Recombinant Human Scavenger receptor cysteine-rich type 1 protein M130 protein (801-940AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Scavenger receptor cysteine-rich type 1 protein M130 (Hemoglobin scavenger receptor) (CD antigen CD163) [Cleaved into: Soluble CD163 (sCD163)], CD163, M130
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	CD163
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.