



RHOT2 Antibody, Biotin conjugated

Product Code	CSB-PA810289LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IXI1
Immunogen	Recombinant Human Mitochondrial Rho GTPase 2 protein (84-300AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mitochondrial Rho GTPase 2 (MIRO-2) (hMiro-2) (EC 3.6.5.-) (Ras homolog gene family member T2), RHOT2, ARHT2 C16orf39
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	RHOT2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.