



FGFBP3 Antibody, FITC conjugated

Product Code	CSB-PA819447OC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TAT2
Immunogen	Peptide sequence from Human Fibroblast growth factor-binding protein 3 protein (192-209AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Fibroblast growth factor-binding protein 3 (FGF-BP3) (FGF-binding protein 3) (FGFBP-3), FGFBP3, C10orf13
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	FGFBP3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.