



GANC Antibody, HRP conjugated

Product Code	CSB-PA837451LB01HU
Abbreviation	GANC
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TET4
Immunogen	Recombinant Human Neutral alpha-glucosidase C protein (171-318AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Relevance	Has alpha-glucosidase activity.
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Neutral alpha-glucosidase C (EC 3.2.1.20), GANC
Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GANC