



TNFAIP8L1 Antibody, HRP conjugated

Product Code	CSB-PA837859LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WVP5
Immunogen	Recombinant Human Tumor necrosis factor alpha-induced protein 8-like protein 1 protein (1-186AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tumor necrosis factor alpha-induced protein 8-like protein 1 (TIPE1) (TNF alpha-induced protein 8-like protein 1) (TNFAIP8-like protein 1) (Oxidative stress-regulated gene-beta) (Oxy-beta), TNFAIP8L1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	TNFAIP8L1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.