



TNFAIP8L1 Antibody, Biotin conjugated

| | |
|----------------------------|---|
| Product Code | CSB-PA837859LD01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q8WVP5 |
| Immunogen | Recombinant Human Tumor necrosis factor alpha-induced protein 8-like protein 1 protein (1-186AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Biotin |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Tumor necrosis factor alpha-induced protein 8-like protein 1 (TIPE1) (TNF alpha-induced protein 8-like protein 1) (TNFAIP8-like protein 1) (Oxidative stress-regulated gene-beta) (Oxy-beta), TNFAIP8L1 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Cell Biology |
| Target Names | TNFAIP8L1 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |