



EDEM1 Antibody, HRP conjugated

Product Code	CSB-PA838793LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92611
Immunogen	Recombinant Human ER degradation-enhancing alpha-mannosidase-like protein 1 protein (46-128AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	ER degradation-enhancing alpha-mannosidase-like protein 1, EDEM1, EDEM KIAA0212
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	EDEM1