



ITCH Antibody, Biotin conjugated

Product Code	CSB-PA839346HD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96J02
Immunogen	Recombinant Human E3 ubiquitin-protein ligase Itchy homolog protein (29-303AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	E3 ubiquitin-protein ligase Itchy homolog (Itch) (EC 2.3.2.26) (Atrophin-1-interacting protein 4) (AIP4) (HECT-type E3 ubiquitin transferase Itchy homolog) (NFE2-associated polypeptide 1) (NAPP1), ITCH
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	ITCH
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.